

I.U.D Procedure Tray

Convenient, disposable, sterile

Catalog # 926/IUD


Contains components for Intra Uterine Device (IUD) insertion procedure.

Pedicat™ Antiseptic gel: lubricates speculum and coats all around vaginal walls, creating safe aseptic environment.

Endofocus Forceps: including gauze ball (Tupfer) for application of antiseptic.

Tenaculum Forceps: For cervical traction. Unique design grasps cervix in a single "snap-close" movement, minimizing pain.

Hysterometer/Uterine Sound: Flexible and foldable round tipped rod with round front edge and moving disk to sound the uterus depth.



Endofocus Clip-Cut: Cuts and removes IUD guiding thread. Unexposed safe blades in "closed" position.


Packaging: 20 trays per case.

Instructions for Use

⚠ The use of this device is limited to those qualified in its use

1. Peel open sealed sterile pack
2. Remove yellow waste disposal bag (can be hung from edge of work surface).
3. Inject/pour *antiseptic solution** into small compartment in the tray.
4. Don *sterile gloves**.
5. Follow instructions printed on the Pedicat™ sachet to lubricate *speculum** with antiseptic lubricating gel.
6. Use speculum in conventional manner by turning it 90° after opening.
7. Grip gauze ball by closing **Endofocus Forceps** jaws, exerting thumb pressure on the rear mechanism (see illustration "A").
8. Using **Endofocus Forceps**, immerse gauze ball in antiseptic solution and disinfect cervical area.
9. Grip cervix edge in one "snap close" movement with **Endofocus Tenaculum Forceps**.
10. While withdrawing **Endofocus Tenaculum Forceps** to align the uterus, sound uterus depth by carefully inserting the **Hysterometer/Uterine Sound**.
11. Adjust I.U.D and **Hysterometer/Uterine Sound** for matching length.
12. Insert I.U.D.
13. To remove **Endofocus Tenaculum Forceps**, open jaws by pulling "trigger" (see illustration "B").
14. Close the **Endofocus Clip-Cut** jaws by maintaining thumb pressure on the rear mechanism (see illustration "C"), insert the device into the vagina, in "closed jaws" position, when close to thread, release thumb pressure, insert thread between jaws and retain pressure to cut, remove **Endofocus Clip-Cut** in closed position with remains of thread.
15. If necessary, excess guiding thread can be removed by using the **Endofocus Forceps**.

* Not included in tray



Pollak Ltd. Hamelacha 20, Rosh Ha'ayin, Afek, 48091 Israel.

Tel: 972 3 9022885 Fax: 972 3 9024531 euroband@hotmail.com

www.euroband.com

EUROBAND

