

SilvaSac™


Central Catheter Waterproof Dressing

Catalog # 1144

Sterile waterproof dressing & pouch, protects exit site of Oncology/Hemodialysis Central Catheters; Holds and preserves catheter tubes and connectors and keeps the exit site and catheter tubes dry*.

Includes:

- **Silvalpad™ – absorbent antibacterial dressing.**
- **Diapouch™ - Soft, elastic waterproof pouch**
(Outer size: 17cm X 8cm) with flap cover for Catheter Connectors.
- **Adhesive tape 10 cm. X 12 cm.**



***Note: The pouch is completely water-resistant. However, pouch flap should completely cover the absorbent antimicrobial dressing with no extruding edges. The adhesive tape should be applied flatly with no kinks and cover the entire flap.**

If well positioned, this pack provides fair protection against water penetration during showering.

Pollak Ltd. Hamelacha 20, Rosh Ha'ayin, Afek, 48091 Israel.
Tel: 972 3 9022885 Fax: 972 3 9024531 euroband@hotmail.com
www.euroband.com

EUROBAND

