



DIABAND™

Prevents blood spurts during
Hemodialysis Procedures.



A soft, pliable pressure armlet, fitted with a highly absorbent non-adherent expanding hemostatic dressing.



Adjustable Ratchet locking system secures a firm, comfortable grip around the arm.

No need for manual pressure following needle/sheath removal.

Pollak Ltd. Hamlacha 20, Rosh Ha'ayin, Afek, 48091 Israel.
Tel: 972 3 9022885 Fax: 972 3 9024531 euroband@hotmail.com
www.euroband.com

EUROBAND

