

# Diapouch™

Catalogue # 01148

## WATERPROOF POUCH HOLDS AND PRESERVES CATHETER TUBES


Sterile, Soft, Elastic Waterproof Pouch  
Outer size: 19X8cm, inner size: 15X7cm.

Protects and manages entry/exit sites of the catheter and its tubing and connectors.

Holds and preserves catheter tubes and connectors of Hemodialysis/Peritoneal/Oncology catheter.

Rear Adhesive strip prevents entry of water during showering.

Provides fair protection against water penetration into the pouch and around the exit site during showering \*



CATHETER IN DIAPOUCH™



CATHETER PRESERVATION WATERPROOF PACK\*

**CATHETER PRESERVATION PACK – DIAPOUCH™** Catalog No. 01148:

**STERILE**, packed in flexible transparent peel-open pouch.

**Contains:**

- **DIAPOUCH™** -Waterproof Pouch, outer size: 19 X 8 cm with additional Rear Adhesive Strip preventing entry of water between skin & pouch.
- **ADHESIVE PLASTER** air permeable, Adhesive in strips, size: 10 X 12 cm.

**\*Note:** *The pouch is completely water resistant. However, pouch flap should completely cover the absorbent dressing around the exit site, with no extruding edges. The adhesive tape should be applied flatly with no kinks and it covers the entire flap. Ensure that the skin is completely dry before adhesion!*

**Pollak Ltd.** Hamlacha 20, Rosh Ha'ayin, Afek, 48091 Israel.

Tel: 972 3 9022885 Fax: 972 3 9024531 euroband@hotmail.com

[www.euroband.com](http://www.euroband.com)

**EUROBAND**

