

POLLAK (INTERNATIONAL) LTD

OFFICES & PLANT: 20 HAMLACHA STREET
NEW INDUSTRIAL ZONE,
ROSH HA'AYIN 48091

EUROBAND
AFEK
ISRAEL VAT No. 510662091

TEL: 972 3 9022885 FAX: 972 3 9024531
e-mail: euroband@hotmail.com
Internet site: www.euroband.com

Formula for Petipad™ Material

	Purpose of Combination	Ingredients	Raw Material Spec.		Amount
Pad	Absorption	Viscose (Rayon)	B.P.C U.S.P E.P	DIN 61210 DIN 53854 DIN 53855	79.97032%
	Strengthens when wet	Polyester ***	U.S.P E.P	DIN 53857 DIN 53923 DIN 53887	19.99258%
	Prevents debris	Polyurethane ****	B.P.C	**	0.03709%
Thread	Textile Thread	Polyester	BSI F 134 : 1989	U.S. Fed. Spec. DD-P-0054a	280 dctx 60%
	The coating of the thread	Barium Sulphate	U.S.P & E.P *****	And B.P.C	39.9999%
		P.V.C *****	ASTM D 1784-92 JIS K 6772		
	Coloring Material	Phtalocyanine Blue *****	U.S F.D.A C.F.R 21&178.3297 B.P.C 1973		0.0001%

- * Thread is the same as the X-Ray Detectable guiding thread raw material specification in U.S. Fed Spec. DD-P-0054a and/or B.P.C 1973 page 627.
- ** D.I.N = German Standards, Covering the materials, production and tests.
- *** Polyester = Same as Polyester Resin specified in Reagent Specifications in U.S.P, and constructed as per DIN 61210
- **** Polyurethane = Same as Polyurethane in The Pharmaceutical CODEX 11th ed. Page 903
- ***** E.P = European Pharmacopoeia.
- ***** P.V.C = Same as Polyvinyl Chloride used in plastics in E.P, or Polyvinyl Chloride (resins) in ASTM, DIN, JIS, AFNOR, KS etc.
- ***** Phtalocyanine blue (Blue Irgalithe BLP produced by Ciba Geigy, Basel Swiss.) Authorized by F.D.A Regulation CFR 21 § 178.3297.

We certify that the above Composition is true.

POLLAK (International) Ltd.

M. Porat
Michael Porat,
Managing Director